

FAST AND EFFICIENT COIL CLEANING



**COILSHINE®-T
DROP-IN
TABLET SYSTEM**



U.S. PATENT NO. 7,841,351

CC-201T – The Answer to Cleaning Coils in Place

The CoilPro® CC-201T, is here to make your on-the-go cleaning of coils a lot quicker and easier. Whether you need to clean coils in mini-split systems, PTAC's, or interior air handlers - CC 201T has got your back.

This system combines a custom pump system, portable design and unique tableted cleaner/bio growth inhibitor in one easy-to-use machine. The CC-201T features an innovative, safe-on-coils, pressurized water system that delivers 200 PSI of cleaning power at 0.9 GPM to remove dirt, dust, and debris from a variety of coils.

It also offers the convenience of a built-in storage tank that holds 8 gallons of water, or you can attach it to a water source for continuous cleaning. Using the unique CoilShine-T tablet system, the CC-201T delivers all the powerful cleaning and protection capabilities of our CoilShine liquid cleaner at a fraction of the size and weight. No more lugging around gallons of cumbersome and bulky coil cleaner.

FEATURES

- Complete, industrial coil cleaning system
- Uses unique CoilShine®-T tablet system (cleaner and bio protectant)
- Powerful, yet safe 200 PSI, 0.9 GPM coil cleaning pressure with 8 gallon tank or continuous water supply
- Distinctive design accommodates all accessories

CC-201T: CC-201T coil cleaner, including spray gun, 25' hose, 18" wand with 3 nozzles. **Call for Price**

CC-201T	
Input Power	Standard: 1.6 Amps, 115 VAC , 50/60 Hz Optional: 0.8 Amps, 230 VAC, 50/60 Hz
Power Supply Cord	21', 14-3 NEMA 515
Water Requirements	On Board or minimum 20 psi @ 1.5 GPM
Output Pressure	130-200 PSI
Output Flow	0.5-0.9 GPM
Pump	Diaphragm
Water Tank Capacity (wet)	8 gallons
Dimensions	44" high x 14.5" wide x 27" length
Weight	65 lbs. dry

**RECOMMENDED FOR
CC-201T & CC-200:
COILSHINE-T**

**RECOMMENDED FOR
CC-200:
COILSHINE-T-BIO**

Find both on pages 59 & 164



Call: 1-800-333-7467
Visit: www.goodway.com
 Goodway Technologies Corporation
 420 West Avenue, Stamford, CT 06902-6384 U.S.A.
 Phone: 203-359-4708 • Fax: 203-359-9601 • goodway@goodway.com